

**REPORT of an EVALUATION**

Microbac

250 West 84<sup>th</sup> Drive

Merillville, IN 46410

IMS LAB #18-B-0137

On

March 29<sup>th</sup>, 2019

By

James B. Lairson

Laboratory Evaluation Officer

Indiana State Board of Animal Health, Dairy Division

1202 East 38<sup>th</sup> Street, Discovery Hall, Suite 100

Indianapolis, IN 46205

Date of Last Evaluation: 03/2017.

Prior Procedures (IMS Code): 2, 5, 20, 21, 22.

Current Procedures (IMS Code): 2, 5, 20, 21, 22.

Prior Lab Status: Fully Accredited.

Present Lab Status: Fully Accredited.

Changes to IMS List: Update Expiration-- Expires 03/2021.

A copy of the Grade "A" Milk Laboratory Evaluation Request and Agreement Form is signed and on file.

**DEVIATIONS AND CORRECTIVE ACTIONS**

**Cultural Procedures- General Requirements 2400 Rev. 4-19.**

<b>Item</b>	<b>Method</b>
<b>3c2.</b>	The thermometers associated with the incubators (16 and 6) were calibrated at 35°C. The above mentioned thermometers need to be calibrated at the temperature of use — In this case 32±1°C.